SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

For

DENTAL IMPLANT SYSTEM (TRADE NAME: RADHEX IMPLANTS)

COMPANY: INMET GARNICK, S.A.

ADDRESS: C./ Del Plástico 5 – Nave 12 Polígono Industrial Miralcampo – Azuqueca De Henares –
(Guadalajara)

General Index of Contents

INTRODUCTION 4

ABREVIATIONS..... 4

THE FOLLOWING INFORMATION IS INTENDED FOR USERS/HEALTHCARE PROFESSIONALS 6

1 The identification of the Device and the Manufacturer, including the basic UDI-DI and, if already issued, the SRN 6

1.1 *Device trade name..... 6*

1.2 *Manufacturer’s name and address 6*

1.3 *Manufacturer’s single registration number (SRN)..... 6*

1.4 *BASIC-UDI-DI..... 6*

1.5 *Medical device nomenclature description/text 6*

1.6 *Class of Device 6*

1.7 *Year when the first certificate (CE) was issued covering the device 7*

1.8 *Authorised representative if applicable, name and the SRN 7*

1.9 *NB’s name (the NB that will validate the SSCP) and the NB’s single identification number 7*

2 The Intended Purpose of the device and any Indications, Contraindications and Target Populations 7

2.1 *Intended Purpose..... 7*

2.2 *Indication(s) and target population(s)..... 7*

2.3 *Contraindications and/or limitations..... 8*

3 A description of the device, including a reference to previous generation(s) or variants if such exist, and a description of the differences, as well as, where relevant, a description of any accessories, other devices and products, which are intended to be used in combination with the device10

3.1 *Description of the device 10*

3.2 *A reference to previous generation(s) or variants if such exist, and a description of the differences. 10*

3.3 *Description of any accessories which are intended to be used in combination with the device. 10*

3.4 *Description of any other devices and products which are intended to be used in combination with the device. 11*

4 Information on any residual risks and any undesirable effects, warnings and precautions.....12

4.1 *Residual risks and undesirable effects 12*

4.2 *Warnings and precautions..... 13*

4.3 *Other relevant aspects of safety, including a summary of any field safety corrective action (FSCA including FSN) if applicable..... 16*

5 The summary of clinical evaluation as referred to in Annex XIV, and relevant information on post-market clinical follow-up.17

5.1 *Summary of clinical data related to equivalent device, if applicable. 17*

5.2 *Summary of clinical data from conducted investigations of the device before the CE marking, if applicable. 17*

5.3 *Summary of clinical data from other sources, if applicable..... 17*

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5.4 *An overall summary of the clinical performance and safety.* 19

5.5 *Ongoing or planned post-market clinical follow-up.* 21

6 **Possible Diagnostic or Therapeutic Alternatives.**21

7 **Suggested profile and training for users.**22

8 **Reference to any harmonized standards and CS applied**.....23

9 **Revision History**.....24

THE FOLLOWING INFORMATION IS A SUMMARY OF THE SAFETY AND PERFORMANCE OF THE DEVICE,25

INTENDED FOR PATIENTS.....25

1 **Device Identification and General Information.**25

Device trade name 25

2 **Intended Use of the Device**.....25

3 **Device description.**28

4 **Risks and warnings.**29

5 **Summary of clinical evaluation and post-market clinical follow-up.**31

6 **Possible diagnostic or therapeutic alternatives**.....32

7 **Suggested training for users.**34

INTRODUCTION

This Summary of Safety and Clinical Performance (SSCP) is intended to provide public access to an up-to-date summary of the main safety and clinical performance aspects related to the device. The SSCP is not intended to replace the Instructions for Use as the primary document to ensure the safe use of the device, nor is it intended to provide diagnostic or therapeutic suggestions to intended users or patients.

This SSCP, which is in two parts, is aimed at:

- Part One: users/health professionals.
- Part Two: patients.

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ABBREVIATIONS

BSI: British Standards Institution

CER: Clinical evaluation report

CIV ID Clinical research identification number, generated by Eudamed for clinical research under the Medical Device Directives

CMR carcinogenic, mutagenic or toxic for reproduction

CS "common specifications" as defined in the MDR EU European Union

C/: Street, address

EMDN: European Medical Device Nomenclature

Eudamed: European database for medical devices

FSCA Field Safety Corrective Action

FSN Field Safety Notice

GMDN: Global Medical Device Nomenclature

IFU Instructions for use

ISO: International Organization for Standardization

MDR Medical Device Regulation

NB notified body

PCI: Internal conical platform

PHE: External hexagonal platform

PHI: Internal hexagonal platform

PMCF Post-market clinical follow-up of

PMS Post-market surveillance of

PSUR Periodic Security Update Report

SRN Unique registration number for an economic operator

SSCP Summary of Safety and Clinical Performance

TD Technical Documentation

Basic UDI-DI Code assigned by the manufacturer to a model or family of medical devices for regulatory purposes.

UDI-DI Unique Device Identification: device identifier

URL Uniform Resource Locator (Internet address)

S.A.: Joint stock company

SLD: Monobloc system; (one piece implant)

SRN: Single registration number

BSI: British Standards Institution

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THE FOLLOWING INFORMATION IS INTENDED FOR USERS/HEALTHCARE PROFESSIONALS

1 The identification of the Device and the Manufacturer, including the basic UDI-DI and, if already issued, the SRN

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1.3 Manufacturer's single registration number (SRN)

The SRN of INMET GARNICK as a manufacturer is:
 ES-MF-000004322
 El SRN de INMET GARNICK como fabricante es:
 ES-MF-000004322

1.4 BASIC-UDI-DI

The BASIC-UDI-DI is:
 For PHI model Dental Implants "84356750PHIA43"
 For PHE model Dental Implants "84356750PHE3T"
 For PCI model Dental Implants "84356750PCI3L"
 For SLD model Dental Implants "84356750SLD4L"
 El BASIC-UDI-DI es:
 For Dental Implants model PHI "84356750PHIA43"
 For Dental Implants model PHE "84356750PHE3T"
 For Dental Implants model PCI "84356750PCI3L"
 For Dental Implants model SLD "84356750SLD4L"

1.5 Medical device nomenclature description/text

The EMDN classifications for the manufactured product correspond to the following:
 P01020101 Implantable Prosthetics and Osteosynthesis devices - Dental Implants.
 Las clasificaciones EMDN que comprenden el producto fabricado corresponde a la siguiente:
 P01020101 Implantable Prosthetics and Osteosynthesis devices- Dental Implants.

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1.7 Year when the first certificate (CE) was issued covering the device

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1.8 Authorised representative if applicable, name and the SRN

Having an authorised representative is not applicable.
No es aplicable el tener un representante autorizado.

1.9 NB's name (the NB that will validate the SSCP) and the NB's single identification number

The name of the NB is: BSI Group the Netherlands B.V. Say Building. John M. Keynesplein 9, 1066 EP Amsterdam, Netherlands. The identification number of the Notified Body is 2797.
El nombre del NB es: BSI Group the Netherlands B.V. Say Building. John M. Keynesplein 9, 1066 EP Amsterdam Netherlands. El número de identificación correspondiente al Organismo Notificado es 2797.

2 The Intended Purpose of the device and any Indications, Contraindications and Target Populations

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2.2 Indication(s) and target population(s)

<p>Indication:</p> <p>It is indicated for clinical conditions of the patient in which:</p> <ul style="list-style-type: none"> - There is a single or multiple absence of teeth in the maxilla and/or mandible. - There is a need for the extraction of teeth that cannot be restored, replaced, or treated by other means. - Other solutions may have an undesirable impact on their healthy teeth or adjacent tissues. <p>Fixed, hybrid or removable prosthetic restoration therapy is diagnosed to re-establish aesthetics and masticatory function, through the application of a surgical protocol in one or two phases, with immediate, early or deferred loading protocols, with sufficient primary stability and adequate occlusal load for the chosen technique. They are indicated for the condition of all types of bone density, with bicortical anchorage being possible in cases of reduced bone density to achieve high initial stability.</p> <p>Target population:</p> <p>The dental implant, as a medical device and considered as a system, involving abutments and abutment placement tools, is intended for a patient population with missing or damaged teeth that need to be replaced to re-establish dental aesthetics and function. These devices can be used in patients, regardless of race, gender, or social condition, as there are no solid reasons that indicate that these factors may affect osseointegration. Warnings regarding contraindications for use are listed in the product's instructions for use.</p>
--

Because dental implants are integrated into the bone by the process of ankylosis, they must be used after the bone structure has completed its growth. Patients must therefore be 18 years of age or older, but the professional must ensure that the patient's bone growth is complete. If this has not occurred, due to the biological condition of the maxillary and mandibular bones where the implant is to be placed, this could affect bone growth, dental eruption and therefore have repercussions on the proper functioning of the product.

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2.3 Contraindications and/or limitations

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General:

- Incomplete bone growth.
- Systemic diseases leading to impaired healing or implant tolerance.
- Uncontrolled diabetes: this can affect the healing process.
- Infectious diseases with risks that compromise the safety of the surgical procedure.
- Allergy to titanium and titanium alloy grade 5 (pre-testing recommended)
- Immunosuppressed or immune system pathologies: in the presence of an autoimmune disease such as lupus or rheumatoid arthritis, the use of dental implants is not recommended until the condition is under control.
- Uncontrolled blood dyscrasias.
- Myocardial infarction (for at least 6 months after the event).
- Thrombo-embolic or cerebral vascular events.
- Endocrinological diseases.
- Sub-acute chronic maxillary osteitis.
- Secondary osteoporosis.
- Osteomalacia.
- Parkinson's. (Relative)

- Psychic alterations.
- Epilepsy.
- Pregnancy.
- Patients with lack of motivation for the implantation.
- Excessive tobacco consumption: this can delay healing, increasing the risk of complications. It is recommended that smokers quit before undergoing the procedure.
- Drug addiction.
- Alcoholism.
- Poor general medical condition of the patient.
- Treatment with anticoagulants (Establish guidelines for expert monitoring of haemostasis prior to surgery).
- Valvular prostheses.
- History of bone radiation. Especially head and neck radiation therapy. (Time assessment and inter-consultation)
- Patient on steroid treatment.
- Prolonged treatment with bisphosphonates can compromise the osseointegration of the implant.
- Chemotherapy.

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3 A description of the device, including a reference to previous generation(s) or variants if such exist, and a description of the differences, as well as, where relevant, a description of any accessories, other devices and products, which are intended to be used in combination with the device

3.1 Description of the device

The RADHEX IMPLANTS system consists of dental implants, related abutments and tools.
 A dental implant is a surgical medical device, which is inserted into the jawbone and acts as a replacement for a missing tooth root. A dental implant consists of two parts: the threaded part which is intended to integrate with the bone, to replace the root of the missing tooth and the upper part, i.e. the dental abutment which is intended to stabilise and retain the prosthetic structure: to fix the tooth crown, the bridge or the denture.

Dental implants are marketed in sterile conditions. Dental implants are sterilised by an electron beam radiation process.
 Table for each product family, with their respective Reference Nomenclatures, and Basic UDI-DI:

Device type or family	Device groups (Reference Nomenclature)	BASIC-UDI-DI
PHI - Internal Hex Platform	PHI350 - PHI450 - PHIA350 - PHIA450 - PHIA+350	84356750PHI43
PHE - External Hexagonal Platform	PHE350 - PHE410 - PHEA410	84356750PHE3T
PCI - Internal Conical Platform	PCI230 - PCI250 - PCI280 - PCI300 - PCI340 - PCI350	84356750PCI3L
SLD - Solid Single Piece Implants	SLD c - SLD cm – SLDb - SLD bm	84356750SLD4L

- Raw material: RADHEX IMPLANTS® Implants, depending on the model, are manufactured in:
- PHI, PHE models in all platforms, made of commercially pure medical grade 4 titanium (ISO standard 5832-2).
 - PHIA models in all platforms, made of medical grade 5 titanium alloy: Ti-6Al-4V (ISO standard 5832-3).
 - PCI platform models, Ø 2.80 mm and 3.50 mm, made of commercially pure grade 4 medical titanium (ISO 5832-2 Standard).
 - Models, PCI platform Ø 2.30 mm, Ø 2.50 mm, Ø 3.00 mm, Ø 3.40 mm
 - Single-body SLD models, made of medical grade 5 titanium alloy: Ti-6Al-4V (ISO standard 5832-3).
 - The healing caps contained in the package of two-piece implants are made of medical grade 5 titanium alloy: Ti-6Al-4V (ISO standard 5832-3).

3.2 A reference to previous generation(s) or variants if such exist, and a description of the differences.

RADHEX IMPLANTS® dental implants do not present any new features compared to similar products on the market, within the current state-of-the-art. The constructions of all these products do not have any major, critical innovations or modifications and are generally accepted at this point in time as state-of-the-art constructions (current level of technology). No clinically relevant changes were made to any of the devices.

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3.3 Description of any accessories which are intended to be used in combination with the device.

Health product accessory: “an article which, without being a medical device itself, is intended by its manufacturer to be used in conjunction with one or more such devices, specifically to enable the medical device(s) to be used in accordance with its intended purpose or to contribute specifically and directly to the medical functionality of the medical device(s) for its intended purpose”.
 Not applicable, as our product does not contain articles that are not themselves medical devices but are intended by their manufacturer to be used in conjunction with our products.

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3.4 Description of any other devices and products which are intended to be used in combination with the device.

The products that are intended to be used in combination with our RADHEX IMPLANTS® dental implants are other complementary products such as abutments and the associated instruments for their proper placement. The indicated products to be used in combination that are not covered by this report are:

- Abutments
- Instrumental

Complementary products may or may not work depending on whether they have the same type of connection or not, for ease of understanding the connection codes used for our products are summarised in the table below:

Types and Correspondence of Dental Implant Connections - Abutment - Screwdrivers			
Implant Code	Corresponding Platform Dimension	Abutment Connection Code	Screwdrivers and Microscrews Tool Connection Code
PHI 350 - PHIA 350 - PHIA+ 350	Two-piece implant: Internal Hexagonal Platform Ø 3.50 mm	IS	HI
PHI 450 - PHIA 450	Two-piece implant: Internal Hexagonal Platform Ø 4.50 mm	IM	
PHE 350	Two-piece implant: External Hexagonal Platform Ø 3.50 mm	ES	ES
PHE 410	Two-piece implant: External Hexagonal Platform Ø 4.10 mm	EM	EM
PCI230	Two-piece implant: Internal Conical Platform Ø 2.30 mm	CS	CS
PCI250	Two-piece implant: Internal Conical Platform Ø 2.50 mm	CS+	CS+
PCI 280	Two-piece implant: Internal Conical Platform Ø 2.80 mm	CM	CM
PCI 300	Two-piece implant: Internal Conical Platform Ø 3.00 mm	CM+	CM+
PCI 350	Two-piece implant: Internal Conical Platform Ø 3.50 mm	CL	CL
PCI340	Two-piece implant: Internal Conical Platform Ø 3.40 mm	CL+	CL+
SLD cm - SLD bm	Solid Single Piece Implant: PMU platform Ø 4.50 mm	PMU	PMU
SLD c - SLD b	Solid Single Piece Implant: Cementable abutment	PMU	PMU

Other medical devices that will be required for the correct placement of dental implants, depending on the technique used, will include, for example, the tools that can be used to prepare the cavity bed where the implant will be placed, by way of example, some of these tools will include the following:

FRL: Lanceolate drill	FRP Pilot Drill	FRF Final Drill	TRJ Threadformer	BIC Circular scalpel
EXP Expander for osteotomy	DIN Torque wrench	PIN Gauge for parallelism, depth and space	PRC Counter angle extender	DTC Contra angle screwdriver

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4 Information on any residual risks and any undesirable effects, warnings and precautions.

4.1 Residual risks and undesirable effects

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Possible risk of patients ingesting or aspirating small devices, so the user should take the appropriate precautionary measures.

UNDESIRABLE SIDE-EFFECTS

Implant surgery or prosthetic intervention may result in:

- Temporary discomfort (pain).
- Bruising and ecchymosis.
- Inflammation of the surgical field. Tissue hyperplasia.
- Localised or systemic infections.
- Fracture of the implant, bone or dental prosthesis.
- Perforation of the maxillary sinus, inferior dental canal and nostrils.
- Accidental introduction of an implant in the maxillary sinus, or nostrils. inflammatory or infectious sinus pathologies. Oroantral or oronasal fistula
- Bone loss due to resorption or recession of soft tissues.
- Deformation or breakage of the implant or abutment due to improper handling.
- Lodging or compression of an implant in the inferior dental canal.

- Unsuccessful implant osseointegration.
- Paraesthesia, dysaesthesia or chronic pain associated with the implant.
- Damage to teeth or neighbouring anatomical structures.
- Intolerance, aesthetic problems, lip bite, bruxism or phonetic alterations.

Immediately after the placement of a dental implant or abutment, activities that require considerable physical exertion should be avoided. Possible complications after dental implants are temporary symptoms: pain, swelling, bleeding, speech difficulties and gingival inflammation.

RADHEX IMPLANTS® implants have passed the ISO 14801 mechanical stress tests. However, they must be properly planned in order to determine the prosthetic loads on them and to reduce the risk of fatigue fracture.

Fracture of implants or their abutments is usually related to loads or forces beyond the strength that the implants or their components can withstand.

Situations that can lead to excessive loading: abutment angles greater than 30 degrees, inadequate number of implants, inadequate implant length or diameter, excessive bridge or crown length, incomplete seating of abutments, occlusal interference, functional problems (e.g. bruxism), improper fit of the prosthesis, lack of passive fit, improper dental laboratory procedures, and injuries caused by habits or accidents suffered by the patient.

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4.2 Warnings and precautions**WARNINGS**

The reuse and/or reprocessing of the single-use implant or medical device creates a potential risk of infection for the patient or user. Contamination can lead to the injury, illness or death of the patient.

The RADHEX IMPLANTS® implant system should not be used in patients or situations presenting contraindications, nor in situations compromising the success of the treatment.

In addition to following the mandatory precautions for any surgical procedure, such as aseptic conditions, the actual lengths of the drills must be identified to determine the working depth with respect to the X-ray measurements, because of the risk of permanent injury to nerves or other vital structures. Drilling deeper than indicated during surgery in the lower jaw may cause permanent sensory disturbances in the lower lip and chin or bleeding in the floor of the mouth. RADHEX IMPLANTS® FRF end drills have 1 mm more extension than the length of the implant to ensure adequate cavity instrumentation, but this must be taken into account when establishing adequate surgical safety margins to avoid injury to vital structures.

Do not use if the package is damaged or opened. Do not use RADHEX dental implants after the expiry date stated on the packaging. Work carefully, maintaining distances and safety margins, avoiding damage to important anatomical structures.

Sterile handling is essential, as contamination can cause infection, and any contact by the medical device with foreign substances prior to use must be avoided.

The product identification label allows each RADHEX IMPLANTS® implant to be tracked, therefore the label should be attached to the traceability card as well as to the patient's medical record.

Instructions for use are provided separately for implants, abutments and instruments or tools, and the appropriate instructions should be followed with respect to the products to be used. For greater understanding and detail you should consult the sections on: Contraindications, warnings, precautions, recommendations, and side effects of RADHEX IMPLANTS®, for abutments and tools.

PRECAUTIONS

Implant surgery and rehabilitation is a complex and highly specialised procedure, and its success cannot be fully guaranteed. For surgical procedure information, see surgical procedure and prosthetic information.

RADHEX IMPLANTS® dental implants, abutments and tools should not be used in locations or situations other than those specifically indicated. RADHEX IMPLANTS® dental implants are intended to be used only with the surgical instruments and abutments belonging to the system, failure to comply with this principle may result in mechanical instrument failure or an unsatisfactory treatment outcome. Gloves must be sterile and powder-free.

Before surgery: Preliminary examination

Patient selection is the responsibility of the professional, who must check that the patient has a suitable case for intraosseous implants and their system components. A thorough physiological and psychological evaluation should be performed, followed by a clinical and radiological examination of the patient prior to the surgical procedure to determine the patient's suitability for treatment.

Any pathological alteration could negatively influence the bone integration or fixation process. Mistakes in treatment planning can lead to implant and abutment failures, especially if the indicated limitations of use or the working procedure are not observed. It is important to determine the appropriate number of implants according to the functional loads, their correct arrangement with respect to the axially of the occlusal forces, avoiding an excessive angle or inclination, their correct distribution and distance between implants, in order to ensure an adequate distribution of the masticatory loads. It is important to take into account the height of the prosthesis, which can represent a significant mechanical risk, especially in single restorations or short restitution sections. Failure to respect these parameters results in the risk of the fracture, deformation or loosening of implants or other prosthetic components.

An adequate anamnesis, a clinical, analytical and radiographic study, must be carried out to confirm that the patient is in a general state of health that does not contraindicate the placement of implants, and the oral cavity must be explored to check its condition and verify that there is no pathology of any kind in the soft or hard tissues.

It is important to check the bone volume or mass and any systemic or local factors that may affect healing. A deficit in hard or soft tissue prior to surgery can lead to an undesirable aesthetic result or unfavourable implant angulations.

Close collaboration between the surgeon, the restorative dentist and the dental laboratory technician is essential at all times for the success of implant treatment. Implant placement and prosthetic design must be tailored to the individual patient's conditions. In case of bruxism, other parafunctional habits or unfavourable intermaxillary relations, a re-evaluation of the treatment option should be considered.

During surgery:

It is recommended that an adequate study of the specific case be carried out, which is the sole responsibility of the acting professional. This is to avoid the risks of prosthetic overload, which must be taken into account, among other things:

- The distribution of possible functional forces as occlusal loads.
- Number of implants.
- Diameter, length and inclination of the implants.
- Transgingival height of prosthetic retainers or retention abutments.
- Possible extension of the arm or cantilever, which may cause a lever effect.

In general, during intraoral use of our products, measures must be taken to prevent aspiration.

During the intervention, the facilities and hygiene infrastructure must be guaranteed. It is necessary to have the appropriate instruments and materials according to the surgeon's needs and indications, depending on the technique to be used and the treatment. The general condition, sterility and especially the sharpness of the instruments must be checked for correct performance.

In the case of implant placement, the surgeon's assessment of the bone quality and primary or initial stability will determine when the implants can be loaded.

Immediate loading is only recommended in the mandible (always splinting implants). Immediate loading is not indicated in non-splinted implants or in cases of low bone density or low implant stability. High torques can cause ischaemia, fracture and compromised bone integration, as well as alterations in the implant connection platform or implant driver. Before attaching the prosthetic component to an implant, the implant must be able to withstand the recommended prosthetic torque. For immediate function, the implant should be able to withstand a torque equivalent to the tightening torque of the inner microscrews. When placing an abutment or cover cap, any hard or soft tissue that may be trapped in the implant connection site must be carefully removed.

It is important to understand that bone can accept surgical trauma, but not thermal trauma; a rise of five degrees Celsius due to friction will lead to bone necrosis and failure of the implant's bone fixation.

After surgery and during prosthetic fitting:

To ensure a good long-term treatment result, it is advisable to carry out a thorough and regular follow-up of the patient after the implant treatment, as well as to inform him/her about the appropriate oral hygiene procedures. Radiotherapy on patients with dental implants should be planned and prescribed with extreme caution by health care professionals to avoid possible complications, and the patient should be adequately informed about the associated risks.

Prosthetic design and fit are essential to the success of implant rehabilitation, and any flaws in this area can lead to treatment failure.

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4.3 Other relevant aspects of safety, including a summary of any field safety corrective action (FSCA including FSN) if applicable

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To date, no other relevant safety aspects have been included nor have any FSCAs and/or FSNs been carried out on our products on the market.

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5 The summary of clinical evaluation as referred to in Annex XIV, and relevant information on post-market clinical follow-up.

5.1 Summary of clinical data related to equivalent device, if applicable.

N/A Not applicable to our product, as no equivalence is claimed.
 N/A No es aplicable a nuestro producto, ya que no se reclama la equivalencia.

5.2 Summary of clinical data from conducted investigations of the device before the CE marking, if applicable.

N/A Not applicable to our product, as no clinical investigations were conducted for the devices prior to CE marking.
 N/A No es aplicable a nuestro producto, ya que no se realizaron investigaciones clínicas para los dispositivos antes del marcado CE.

5.3 Summary of clinical data from other sources, if applicable.

Data collection for the PMCF is carried out on post-market data. The annual study was initiated with the aim of validating our products on the market, obtaining data on the success rate of the products during day-to-day practice, the current state-of-the-art, their usability and intended use. All devices sold were included in the study without exception.

The percentage of successful-osseointegrated products per year is shown below:

	2018	2019	2020	2021	2022
Total	99.4121%	99.1579%	98.7110%	98.6625%	98,492%

During the years of follow-up, no new adverse events or side effects have been directly related to our product. Therefore, those indicated in the instructions for use are maintained.

No alerts or withdrawals from the market have been registered so far.

Clinical data on the devices in question were collected through a Post Market Clinical Follow-Up clinical experience study for the CE mark for the INMET GARNICK implant system. The study plan was initiated with the aim of validating the INMET GARNICK dental implant system for its usual practice for products in relation to the state-of-the-art and its proposed intended use. All devices belonging to INMET GARNICK were included in the study without exception.

The primary aim of the study was to obtain the survival rate of the implants 5 years after implantation. The secondary objectives are to determine, evaluate and quantify possible complications related to the implantation procedure in the short, medium and long term.

In order to avoid potential manufacturer bias as much as possible, no patient inclusion or exclusion criteria were applied. These are the criteria that are usually encountered in the day-to-day practice of professionals.

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All the implants will undergo a 5-year follow-up. The data collected so far showed excellent clinical results: a cumulative survival rate for INMET GARNICK dental implants according to the platform of:

- PHI: 97.23% (5 years after placement)
- PHE: 97.89% (5 years after placement)
- PCI: 96.55% (3 years after placement)
- SLD: 97.89% (2 years after placement)

The complications presented had the following results: pain (0-1.44%), tenderness (0-1.63%), swelling (3.06-19.42%), bleeding (0.51-1.63%), suppuration (0%), bone loss from 0mm to 1mm (7.65-19.51%), bone loss from 1mm to 2mm (4.32-6.12%), bone loss 2mm to 3mm (0.41-2.16%), bone loss 3mm to 4mm (0-1.44%) fractured (0-1.63%), connection incidents (0-3.60%).

No adverse events, side effects or unforeseen or unexpected complications were reported during the evaluation.

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5.4 An overall summary of the clinical performance and safety.

RADHEX IMPLANTS® has used clinical experience with the product as a source of direct clinical data on its manufactured devices, and these are monitored as part of the PMCF.

This makes it possible to maintain constant monitoring of the product's clinical performance data, and to be sensitive to compliance with the necessary safety and performance requirements for the devices.

We then proceed to the compilation and collation of the data collected.

INMET GARNICK began its clinical record in 2017 as part of its PMCF. During the PMCF study, the safety and performance of the device was closely monitored. The clinical data are inconclusive, as we anticipate data after 5 years of implantation, however, we currently have data between the third and fifth year of implantation for implants depending on the implant connection platform (PHE 5th year of evaluation, PHI 5th year of evaluation, PCI 3rd year of evaluation and SLD 2nd year of evaluation).

Overall, a total of 740 implants were placed (Divided by platform into 142-PHE; 253-PHI; 203-PCI and 142-SLD) in a total of 228 patients (Divided by platform into 54-PHE; 50-PHI; 74-PCI and 50-SLD). The implants were fitted in the period between 2015 and 2020.

The results obtained by the PMCF study so far are in line with the available clinical data pertaining to this group of devices on the market. The results obtained for the survival rate of the implants placed are as follows:

PHI: 97.23% (5 years after placement)

PHE: 97.89% (5 years after placement)

PCI: 96.55% (3 years after placement)

SLD: 97.89% (2 years after placement)

The performance of the medical device is described as the ability to achieve the intended use claimed by the manufacturer.

Therefore, the value of the survival rate of dental implants can be considered as the main indicator of implant performance.

Another indicator of performance is the absence of mechanical failure, such as fracture due to the fatigue of implants or attachments. Consistent with the PMCF study, there were no adverse events, including failure due to fatigue of either the implant or the abutment.

Medical device safety is described as the acceptability of risks versus benefits, when the medical device is used according to the manufacturer's instructions.

Safety concerns regarding dental implant safety are mostly related to complication rates (adverse effects) and risks associated with implant loading procedures. The non-occurrence of severe side effects is an indicator for the clinical safety of implants. As the PMCF study showed no severe side effects for dental implants, it can be considered to be of the highest clinical safety in this aspect during the years of follow-up.

For the PHI platform after 5 years of implantation (253 implants, 50 patients)

The only adverse effects recorded after treatment with the dental implant system were temporary and not severe, including: pain (0.41%), tenderness (1.63%), swelling (3.66%), bleeding (1.63%), suppuration (0%), bone loss of 0 mm to 1 mm (19.51%), bone loss 1 mm to 2 mm (4.88%), bone loss 2 mm to 3 mm (0.41%), fractures (1.63%), connection incidents (1.22%). The significantly low suppuration ratio (0%) confirms that the implant surface is of the highest quality and cleanliness, providing a safe and desirable connection between the bone and the implant and ensuring osseointegration.

For the PHE platform after 5 years of implantation (142 implants, 54 patients)

The only adverse effects recorded after treatment with the dental implant system were temporary and not severe, including: pain (1.44%), tenderness (0%), swelling (19.42%), bleeding (0.72%), suppuration (0%), bone loss of 0mm to 1mm (17.27%), bone loss 1mm to 2mm (4.32%), bone loss 2mm to 3mm (1.44%), bone loss 3mm to 4mm (1.44%), fractures (0%), connection incidents (3.60%). The significantly low suppuration ratio (0%) confirms that the implant surface is of the highest quality and cleanliness, providing a safe and desirable connection between the bone and the implant and ensuring osseointegration.

For the PCI platform after 3 years of implantation (203 implants, 74 patients)

The only adverse effects recorded after treatment with the dental implant system were temporary and not severe, including: pain (0.51%), tenderness (0.51%), swelling (3.06%), bleeding (0.51%), suppuration (0%), bone loss of 0mm to 1mm (7.65%), bone loss 1mm to 2mm (6.12%), bone loss 2mm to 3mm (0.51%), fractures (0%), connection incidents (0%). The significantly low suppuration ratio (0%) confirms that the implant surface is of the highest quality and cleanliness, providing a safe and desirable connection between the bone and the implant and ensuring osseointegration.

For the SLD platform after 2 years of implantation (142 implants, 50 patients)

The only adverse effects recorded after treatment with the dental implant system were temporary and not severe, including: pain (0%), tenderness (0%), swelling (11.51%), bleeding (1.44%), suppuration (0%), bone loss of 0mm to 1mm (15.65%), bone loss 1mm to 2mm (5.76%), bone loss 2mm to 3mm (2.16%), fractures (1.63%), connection incidents (0%). The significantly low suppuration ratio (0%)

confirms that the implant surface is of the highest quality and cleanliness, providing a safe and desirable connection between the bone and the implant and ensuring osseointegration.

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5.5 Ongoing or planned post-market clinical follow-up.

The current PMCF continues to collect information on possible incidents related to the product, obtaining information in the medium to short term, as it is in this period when there is a greatest probability that osseointegration will not occur.

The extended objective of the study is to achieve a 90% osseointegration rate for the implant in the first months after placement. The following questions will be answered with more clinical data:

- Will the target of 90% osseointegration success be achieved?

In the future, the aim is to find out how our products will perform in the medium term with the immediate loading technique applied. The aim of the study is to find out whether the implants inserted using this technique has a survival rate of at least 90% after a period of 5 years.

The following questions will be answered with more clinical data:

- Will a 90% survival rate be achieved through the application of the immediate loading technique?

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6 Possible Diagnostic or Therapeutic Alternatives.

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- Crowns/bridges: this is the restorative treatment of a tooth-supported fixed prosthesis, consisting of replacing the missing teeth, taking support from neighbouring teeth to the edentulous section.

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Advantages: it is a treatment that allows functional and aesthetic recovery, maintaining part of the original dentition.
Disadvantages: It requires healthy teeth close to the edentulous area. In addition, it requires the coronal carving approach of the neighbouring teeth to anchor the prosthetic structure.

- Anchorless prosthesis. This is the restorative treatment of removable prostheses.

Advantages: It is a classic non-surgical treatment for the restoration of teeth.

Disadvantages: Removable prostheses have lower functional efficacy, may cause chewing insufficiency, phonatory insufficiency, soft tissue friction and lack of confidence in patients who are often restricted to soft, easy to chew foods and are not protected from bone loss. It can also cause discomfort due to food debris that gets between the prosthesis and the soft tissues.

- Alternative intraosseous implants. There are many different types of dental implants on the European market.

Advantages: Many of them have proven effectiveness and have been on the market for several years.

Disadvantages: Not all of them are marketed in Spain or offer advantages from a technical point of view.

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7 Suggested profile and training for users.

RADHEX IMPLANTS® implants, abutments and tools may only be used in the dental health care environment by certified licensed professionals who are properly trained to work with them. The intended user is the professional who, due to his or her medical training and legal authorisation, is able to perform surgery in the oral and maxillary cavity. These training courses correspond to qualifications such as dentistry, dental surgeon, dentist, oral surgeon or maxillofacial surgeon.

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8 Reference to any harmonized standards and CS applied.

- ISO 14644-1: Clean rooms and controlled adjoining premises Part 1: Classification of air cleanliness by particle concentration
- EN 1642: Dentistry - Medical devices for dentistry - Dental implants
- MEDDEV 2.7/1 Rev. 4. Guidelines for medical devices. Clinical evaluation: guidance for manufacturers and notified bodies.
- ISO 13485 Medical devices. Quality management system: requirements for regulatory proposals
- ISO 14971 Medical devices - Application of risk management to medical devices
- ISO 14801 Fatigue test for endosseous dental implants
- ISO 10993-1 Biological evaluation of medical devices. Evaluation and testing
- ISO 10993-5 Biological evaluation of medical devices. Tests for in vitro cytotoxicity
- ISO 10993-10 Biological evaluation of medical devices. Skin irritation and sensitisation test.
- ISO 11607-1: Packaging for terminally sterilised medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems
- ISO 11607-2: Packaging for terminally sterilised medical devices - Part 2: Validation requirements for form, seal and assembly processes
- MEDDEV 2.12-1: Guidelines for a monitoring system for medical devices
- ISO 15223-1 Medical devices - Symbols to be used with medical device labels, labelling and information to be provided. Part 1: general requirements
- ISTA 3A: Test procedure packaging performance simulation for 70 kg or less shipping system
- ISO 20417 Medical devices. Information to be provided by the manufacturer.
- ISO 11737-1 Sterilisation of health products. Microbiological methods. Part 1: Determination of the population of micro-organisms in products.
- ISO 11737-2: Sterilisation of health products. Microbiological methods. Part 2: Sterility tests carried out for the definition, validation and maintenance of a sterilisation process.
- ISO 11137-1: Sterilisation of products for health care. Radiation. Part 1: Requirements for the development, validation and routine control of a sterilisation process for medical devices.
- ISO 11137-2: Sterilisation of products for health care. Radiation. Part 2: Establishing the sterilisation dose.
- ISO 11137-3: Sterilisation of products for health care. Radiation. Part 3: Recommendations on the dosimetric aspects of development, validation and routine monitoring.
- IEC 62366-1 Medical devices. Part 1: Application of usability engineering to medical devices.
- ISO 10451: Dentistry. Content of the technical file for dental implant systems.
- ISO 7405: Dentistry. Evaluation of the biocompatibility of medical devices used in dentistry.
- EN 868-2 Packaging materials for sterilised medical devices in their final phase. Part 2: Sterilisation wrap. Requirements and test methods.
- EN 17141: Clean rooms and associated controlled environments. Control of biocontamination.
- ISO 5832-2: Surgical implants. Metallic materials. Part 2: Unalloyed titanium.
- ISO 5832-3: Surgical implants. Metallic materials. Part 3: Wrought alloy from titanium 6-aluminium 4-vanadium
- ISO 14644-1: Salas limpias y locales anexos controlados Parte 1: Clasificación de la limpieza del aire mediante la concentración de partículas
- EN 1642: Odontología Productos sanitarios para odontología Implantes dentales
- MEDDEV 2.7/1 Rev. 4. Directrices para dispositivos médicos. Evaluación clínica: guía para fabricantes y organismos notificados.
- ISO 13485 Productos sanitarios. Sistema de gestión de la calidad: requisitos para propuestas reglamentarias
- ISO 14971 Productos sanitarios - Aplicación de la gestión de riesgos a los productos sanitarios
- ISO 14801 Ensayo de fatiga para implantes dentales endoóseos
- ISO 10993-1 Evaluación biológica de dispositivos médicos. Evaluación y prueba
- ISO 10993-5 Evaluación biológica de productos sanitarios. Prueba de citotoxicidad in vitro
- ISO 10993-10 Evaluación biológica de dispositivos médicos. Prueba de irritación y sensibilización cutánea.
- ISO 11607-1: Envasado para productos sanitarios esterilizados terminalmente Parte 1: Requisitos para los materiales, los sistemas de barrera estéril y sistemas de envasado
- ISO 11607-2: Envasado para productos sanitarios esterilizados terminalmente Parte 2: Requisitos de validación para procesos de conformación, sellado y ensamblado
- MEDDEV 2.12-1: Directrices para sistema de vigilancia en productos sanitarios

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9 Revision History.

SSCP revisión number	Date issued	Change description	Revision validated by the Notified Body	
00	07.10.2021	FIRST ISSUE: Summary of Safety and Clinical Performance of Medical Devices (SSCP)		Yes Validation language:
			X	No (only applicable for class IIa or some IIb implantable devices (MDR, Article 52(4) (2nd paragraph) for which the SSCP is not yet validated by the NB)
01	30.06.2022	Updated 2021 data		Yes Validation language:
			X	No
02	08.09.2022	Review of point 8 from NC 2241147-202209-N7		Yes Validation language:
			X	No
03	28.09.2022	Updated documentation		Yes Validation language:
			X	No
04	17.11.2022	Updated documentation		Yes Validation language:
			X	No
05	16.03.2023	Updated documentation	X	Yes Validation language: English
				No
06	16.06.2023	Updated documentation	X	Yes Validation language: English
				No

INMET GARNICK S.A.	
SSCP	I-1-SSCP-Implants Report
SUMMARY OF SAFETY AND CLINICAL PERFORMANCE	Edition 06

THE FOLLOWING INFORMATION IS A SUMMARY OF THE SAFETY AND PERFORMANCE OF THE DEVICE, INTENDED FOR PATIENTS

Document revision: 04
Date issued: 2023.06.16

This Summary of Safety and Clinical Performance (SSCP) is intended to provide public access to an up-to-date summary of the main safety and clinical performance aspects of the device. The following information is intended for patients or non-professionals. A more extensive summary of its safety and clinical performance prepared for health professionals can be found in the first part of this document.

The SSCP is not intended to provide general advice on the treatment of a medical condition. Contact your health care professional if you have questions about your medical condition or the use of the device in your situation. This SSCP is not intended to replace an implant card or the Instructions for Use to provide information on the safe use of the device.

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1 Device Identification and General Information.

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Basic UDI-DI

PHI "84356750PHIA43"
PHE "84356750PHE3T"
PCI "84356750PCI3L"
SLD "84356750SLD4L"
PHI "84356750PHIA43"
PHE "84356750PHE3T"
PCI "84356750PCI3L"
SLD "84356750SLD4L"

Year when the device was first CE-marked

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2 Intended Use of the Device.

INMET GARNICK S.A.	
SSCP	I-1-SSCP-Implants Report
SUMMARY OF SAFETY AND CLINICAL PERFORMANCE	Edition 06

Intended purpose

This device is intended to replace missing or damaged teeth, by means of fixed, hybrid or removable, single or multiple, partial or total prosthetic rehabilitation techniques, thus restoring dental aesthetics and functionality by anchoring, retaining and stabilising dental prostheses

The dental implant system consists of a main piece, which is the implant itself, and a series of accessory pieces, called abutments, which are connected to it and behave as a functional unit, and the instruments that allow these pieces to be placed and connected to each other, so that they can fulfil the described intended use.

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Indications and target population

Indication:

It is indicated for clinical conditions of the patient in which:

- There is a single or multiple absence of teeth in the maxilla and/or mandible.
- There is a need for the extraction of teeth that cannot be restored, replaced, or treated by other means.
- Other solutions may have an undesirable impact on their healthy teeth or adjacent tissues.

Fixed, hybrid or removable prosthetic restoration therapy is diagnosed to re-establish aesthetics and masticatory function, through the application of a surgical protocol in one or two phases, with immediate, early or deferred loading protocols, with sufficient primary stability and adequate occlusal load for the chosen technique. They are indicated for the condition of all types of bone density, with bicortical anchorage being possible in cases of reduced bone density to achieve high initial stability.

Target population:

The dental implant, as a medical device and considered as a system, involving abutments and abutment placement tools, is intended for a patient population with missing or damaged teeth that need to be replaced to re-establish dental aesthetics and function.

These devices can be used in patients, regardless of race, gender, or social condition, as there are no solid reasons that indicate that these factors may affect osseointegration. Warnings regarding contraindications for use are listed in the product's instructions for use. Because dental implants are integrated into the bone by the process of ankylosis, they must be used after the bone structure has completed its growth. Patients must therefore be 18 years of age or older, but the professional must ensure that the patient's bone growth is complete. If this has not occurred, due to the biological condition of the maxillary and mandibular bones where the implant is to be placed, this could affect bone growth, dental eruption and therefore have repercussions on the proper functioning of the product.

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Contraindications

In general, the insertion of Rhadex System® implants is contraindicated in patients:

- Who are not medically fit to undergo an oral surgical procedure.
- With insufficient bone volume, unless an augmentation procedure can be considered.
- Who cannot be fitted with the appropriate size, number or desirable position of implants to safely support functional or parafunctional loads in the long term.
- Who are allergic or hypersensitive to commercially pure grade 4 titanium or Ti-6Al-4V titanium alloy (titanium, aluminium and vanadium)
- Contraindicated in children, paediatric patients and adolescents who have not completed bone growth.

In detail we mention the following contraindications:

Local:

- Local bone deficiency or atrophic residual bone that makes implant placement difficult or impossible.
- Excessive occlusal parafunction.
- Patients with poor oral hygiene.
- Active or advanced local pathologies: cystic, infectious, tumoural or periodontal, requiring prior treatment
- Paranasal sinus pathologies.

General:

- Incomplete bone growth.
- Systemic diseases leading to impaired healing or implant tolerance.
- Uncontrolled diabetes: this can affect the healing process.
- Infectious diseases with risks that compromise the safety of the surgical procedure.
- Allergy to titanium and titanium alloy grade 5 (pre-testing recommended)
- Immunosuppressed or immune system pathologies: in the presence of an autoimmune disease such as lupus or rheumatoid arthritis, the use of dental implants is not recommended until the condition is under control.
- Uncontrolled blood dyscrasias.
- Myocardial infarction (for at least 6 months after the event).
- Thrombo-embolic or cerebral vascular events.
- Endocrinological diseases.
- Sub-acute chronic maxillary osteitis.
- Secondary osteoporosis.
- Osteomalacia.
- Parkinson's. (Relative)
- Psychic alterations.
- Epilepsy.
- Pregnancy.
- Patients with lack of motivation for the implantation.
- Excessive tobacco consumption: this can delay healing, increasing the risk of complications. It is recommended that smokers quit before undergoing the procedure.
- Drug addiction.
- Alcoholism.
- Poor general medical condition of the patient.
- Treatment with anticoagulants (Establish guidelines for expert monitoring of haemostasis prior to surgery).
- Valvular prostheses.
- History of bone radiation. Especially head and neck radiation therapy. (Time assessment and inter-consultation)
- Patient on steroid treatment.
- Prolonged treatment with bisphosphonates can compromise the osseointegration of the implant.
- Chemotherapy.

The professional is responsible for informing the patient about all possible contraindications, side effects and precautions, as well as for providing maintenance and expert advice if there is any alteration in the results of the implant (e.g. loosening of the prosthesis, inflammation, infection or oozing around the implant, pain or any other unexpected symptom).

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3 Device description.

Device description and material/substances in contact with patient tissues

RADHEX IMPLANTS consists of dental implants, related abutments and tools. A dental implant is a surgical accessory that is placed in the jawbone and acts as a replacement for the root of a missing tooth. A dental implant consists of two parts: the threaded part which is intended to fuse with the bone, to replace the root of the missing tooth and the upper part, i.e. the dental abutment which is intended to fix the tooth crown, the bridge or the denture.

RADHEX IMPLANTS® Implants, depending on the model, are manufactured in:

- PHI, PHE models in all platforms, made of commercially pure medical grade 4 titanium (ISO standard 5832-2).
- PHIA models in all platforms, made of medical grade 5 titanium alloy: Ti-6Al-4V (ISO standard 5832-3).
- PCI platform models, Ø 2.80 mm and 3.50 mm, made of commercially pure grade 4 medical titanium (ISO 5832-2 Standard).
- Models, PCI platform Ø 2.30 mm, Ø 2.50 mm, Ø 3.00 mm, Ø 3.40 mm
- Single-body SLD models, made of medical grade 5 titanium alloy: Ti-6Al-4V (ISO standard 5832-3).
- The healing caps contained in the package of two-piece implants are made of medical grade 5 titanium alloy: Ti-6Al-4V (ISO standard 5832-3).

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Information about medicinal substances in the device, if any

The medical device does not contain medicinal substances.

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Description of how the device is achieving its intended mode of action

Action mode:

Fixation to the bone and replacement of the tooth root. To achieve the intended use, it will be necessary to complement it with an abutment, which will work by attaching the required prostheses to the implant.

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Description of accessories, if any

Does not apply. RADHEX IMPLANTS® does not have accessories for its products.

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4 Risks and warnings.

Contact your healthcare professional if you think you are experiencing side effects related to the device or its use or if you are concerned about risks. This document is not intended to replace a consultation with your health care professional if necessary.

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How potential risks have been controlled or managed

The potential risks of the product have been managed or controlled starting from the design of the product, the production process and the sterilisation of the product, thus allowing the product to arrive in perfect condition until the moment of implantation in the patient and reducing potential risks with the product to a minimum.

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Remaining risks and undesirable effects

In the field of medical technologies, absolute safety can never be guaranteed, and this must be taken into account by the user when deciding on treatment.

Failure to comply with the instructions for use, limitations of use, inappropriate use of the products, failure to take into account warnings, precautions and recommendations, reuse of the device, lack of collaboration between the surgeon, the restorative dentist and the laboratory technician, mechanical failure as the appropriate insertion torques not respected, implantation in contraindicated conditions, can lead to treatment failures and severe damage that can compromise the patient's well-being, health and life.

It should be borne in mind that implant treatment can have negative consequences that carry risks in biological, structural, functional, psychological and mechanical terms, which include the following:

Bone loss, injuries to adjacent anatomical structures, nerve or sensory injuries ranging from numbness to neuralgia, vascular injuries or haemorrhages, local potentially systemic infections, aesthetic and phonetic alterations, failure of integration, prosthetic mobility, platform deformation, implant fracture.

Temporary discomfort after invasive treatment, such as side effects, is common.

The risk of patients ingesting or aspirating small devices is possible and the user should take the appropriate precautionary measures.

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Warnings and precautions

WARNINGS
 The reuse and/or reprocessing of the single-use implant or medical device creates a potential risk of infection for the patient or user. Contamination can lead to the injury, illness or death of the patient.
 The RADHEX IMPLANTS® implant system should not be used in patients or situations presenting contraindications, nor in situations compromising the success of the treatment.
 In addition to following the mandatory precautions for any surgical procedure, such as aseptic conditions, the actual lengths of the drills must be identified to determine the working depth with respect to the X-ray measurements, because of the risk of permanent injury to nerves or other vital structures. Drilling deeper than indicated during surgery in the lower jaw may cause permanent sensory disturbances in the lower lip and chin or bleeding in the floor of the mouth. RADHEX IMPLANTS® FRF end drills have 1 mm more extension than the length of the implant to ensure adequate cavity instrumentation, but this must be taken into account when establishing adequate surgical safety margins to avoid injury to vital structures.
 Do not use if the package is damaged or opened. Do not use RADHEX dental implants after the expiry date stated on the packaging. Work carefully, maintaining distances and safety margins, avoiding damage to important anatomical structures.
 Sterile handling is essential, as contamination can cause infection, and any contact by the medical device with foreign substances prior to use must be avoided.
 The product identification label allows each RADHEX IMPLANTS® implant to be tracked, therefore the label should be attached to the traceability card as well as to the patient's medical record.
 Instructions for use are provided separately for implants, abutments and instruments or tools, and the appropriate instructions should be followed with respect to the products to be used. For greater understanding and detail you should consult the sections on: Contraindications, warnings, precautions, recommendations, and side effects of RADHEX IMPLANTS®, for abutments and tools.

PRECAUTIONS
 Implant surgery and rehabilitation is a complex and highly specialised procedure, and its success cannot be fully guaranteed. Specialised training in established techniques or oral implantology is required. Courses on implants and seminars are highly recommended. For surgical procedure information, see surgical procedure and prosthetic information.
 RADHEX IMPLANTS® dental implants, abutments and tools should not be used in locations or situations other than those specifically indicated. RADHEX IMPLANTS® dental implants are intended to be used only with the surgical instruments and abutments belonging to the system, failure to comply with this principle may result in mechanical instrument failure or an unsatisfactory treatment outcome. Gloves must be sterile and powder-free.

After surgery and during prosthetic fitting:
 To ensure a good long-term treatment result, it is advisable to carry out a thorough and regular follow-up of the patient after the implant treatment, as well as to inform him/her about the appropriate oral hygiene procedures. Radiotherapy on patients with dental implants should be planned and prescribed with extreme caution by health care professionals to avoid possible complications, and the patient should be adequately informed about the associated risks.
 The design and fit of the prosthesis are essential to the success of implant rehabilitation, and any problems in this area can lead to treatment failure.

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Summary of any field safety corrective action, (FSCA including FSN) if applicable

No field safety corrections or field safety notifications have been necessary so far.

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5 Summary of clinical evaluation and post-market clinical follow-up.

Clinical background of the device

Dental implants were introduced in the 1960s and have been used ever since as a treatment option in dentistry. Despite the great variety of different sizes (e.g. diameter, length, shape), thread design and/or surface topography, the mode of action of dental implants (osseointegration, i.e. the ability to fuse with the bone), ensuring fixation in the jawbone, and the intended purpose has remained almost identical over the years. A similar assumption can be made for the surgical implantation procedure, therefore, dental implants can be considered the state-of-the-art treatment for patients with missing teeth. More than 50 years of proven clinical safety and performance related to dental implants and related abutments are widely available in the scientific literature.

RADHEX IMPLANTS has been marketing its products since 2009, the year in which it obtained its corresponding European CE certification. It has 14 years of experience in the market and uses the same materials that have been tested for safety through biocompatibility testing and similar designs as other dental implant systems available on the market.

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The clinical evidence for the CE-marking

RADHEX IMPLANTS has well-established safety and clinical performance, which is based on the biocompatibility of the materials used, established and validated cleanliness and shelf-life stability, and real-life clinical data, based on implant survival rate.

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Safety

Dental implant treatment has a clearly positive risk-benefit profile. The manufacturer has considered the totality of the evidence regarding the extent of the likely benefits and the extent of the likely risks of a device in the risk-benefit information. All residual risks were mitigated to the extent possible until additional risk controls are no longer practicable. All identified hazardous situations have been assessed and all risks have been reduced to an acceptable level based on a risk-benefit analysis. Every residual risk is acceptable, the benefits far outweigh the risks. No additional side effects from residual risks were identified. On the basis of the benefits and risks, it can be concluded that the multiple benefits far outweigh the residual risks. No additional risk control measures are applied and the medical device, i.e. RADHEX IMPLANTS can be considered safe for its intended use.

The manufacturer continuously collects clinical data and information on the safety and clinical performance of the device by conducting a post-market clinical follow-up study. Careful observation of implant survival and short- and long-term side effects is carried out during the study.

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6 Possible diagnostic or therapeutic alternatives.

When considering alternative treatments, it is recommended that you contact your health care professional, who can take into account your individual situation when determining the therapeutic alternative of choice for your oral problem.

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General description of therapeutic alternatives

Dental implants are the last resort on a dentist's spectrum of possibilities for restoring the ability to chew and are commonly used when tooth preservation is no longer a safe and realistic option.

Since abutments do not have independent indications and are intended for use in combination with dental implants already in place, with which they form a functional unit to fulfil the intended purpose, the indication of the implants must be taken into account when considering therapeutic alternatives.

When a dentition deficiency affects the ability to chew, there are several therapeutic alternatives to implants, depending on the severity of the dental defects:

- Conservative treatment. These are traditional dental treatments aimed at preserving the teeth, and are the main alternative as long as it is possible to preserve the tooth and restore its health and functionality through these techniques.

Advantages: Conservative treatments can be applied, maintaining the natural tooth or most of its structure.

Disadvantages: The application of treatments to maintain teeth is not possible when the tooth has been lost, or when the technique does not make it possible to maintain the health and function of the tooth.

- Dental transplant / Autotransplant. The autotransplantation of dental pieces is an alternative technique that consists of using a tooth that has been extracted, to position it in the edentulous area, in a socket created to stabilise it after healing.

Advantages: The resource is obtained from the patient directly.

Disadvantages: A tooth is required for the recipient site, and it is a non-standardised technique, requiring more time and high operator skill which makes it difficult to apply.

- Crowns/bridges: this is the restorative treatment of a tooth-supported fixed prosthesis, consisting of replacing the missing teeth, taking support from neighbouring teeth to the edentulous section.
 Advantages: it is a treatment that allows functional and aesthetic recovery, maintaining part of the original dentition.
 Disadvantages: It requires healthy teeth close to the edentulous area. In addition, it requires the coronal carving approach of the neighbouring teeth to anchor the prosthetic structure.
- Anchorless prosthesis. This is the restorative treatment of removable prostheses.
 Advantages: It is a classic non-surgical treatment for the restoration of teeth.
 Disadvantages: Removable prostheses have lower functional efficacy, may cause chewing insufficiency, phonatory insufficiency, soft tissue friction and lack of confidence in patients who are often restricted to soft, easy to chew foods and are not protected from bone loss. It can also cause discomfort due to food debris that gets between the prosthesis and the soft tissues.
- Alternative intraosseous implants. There are many different types of dental implants on the European market.
 Advantages: Many of them have proven effectiveness and have been on the market for several years.
 Disadvantages: Not all of them are marketed in Spain or offer advantages from a technical point of view.

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The medical device manufactured by Inmet Garnick S.A. comprises surgical anchor fixation systems, with the purpose of providing support, anchoring and stabilisation in the surgical field, bearing the trademark Radhex Implants.